


changes when the phase transition is

Pro diffractometer and Anton Paar

sample holder were compensated using a

to trigonal(II) phase at $T(\text{Lit})=129\pm1^\circ\text{C}$
 ms at 120°C before conversion back to



se formation during cooling, starting with
 ase(II).

green curve shows the transition of
 (III) for a second cooling cycle.

ing cycle is $T(\text{meas})=128^\circ\text{C}$, which is
 ng down the sample (Fig.2), a sharp
 at 120°C . This transition is very
) proceeds slowly and is not complete.
 e at which the transition (I) \Rightarrow (II) takes

Fig. 4

Transition of BaCO_3 from phase(I) \Rightarrow phase
 samples (s#1,s#2) for repeated heating c
 (m#1-3).

The Graph shows the decrease of the phase
 intensity at $\approx 23.4^\circ 2\theta$

The transition temperature of $T(\text{meas})$
 tabulated value.

For subsequent heating cycles the t

The reason for the observed shift in tr
 BaCO_3 powder during the first heating

Because of the shown change in trans
 and its destructive effect on platinum h
XRD temperature calibration.

Temperature calibration w

As an example for a calibration with pl
 shows the calibration curve obtained f




Fig. 6

Calibration curve with phase transitions
 for an XRK900 heating attachment.

Atmosphere in the chamber: air, 1 bar

Calibration of the temperature measur