



# BRITISH CRYSTALLOGRAPHIC ASSOCIATION

Charity Registration Number 284718

## UK Instrument Sensitivity Round Robin

D. J . Taylor. BCA Industrial Group

### Introduction.

The seeds of the Round Robin exercise were sown at the Spring meeting at Cardiff in 1995 during discussions with Ron Jenkins of ICDD. It stemmed from a need to evaluate an individual instrument's performance against other instruments, of similar configuration, and set a sensitivity benchmark against which any instrument could be judged. A similar exercise had been completed in the USA and plans were well advanced for other countries to run the same exercise. At the Industrial Group Committee meeting on 2nd May 1995 it was decided to run the USA exercise in the UK. Ron Jenkins supplied details of the experiments, a Quattro Pro spreadsheet to evaluate the measured data and details of the USA results. The event was advertised in Crystallography News and IG Newsletter and the first registration was received in July 1995. In total 35 instruments were registered for the exercise and 33 tests were completed by January 1997 and are included in this report.


### Basis of Tests.

To ensure that the participant did not introduce any sample preparation error into the tests the NIST SRM1976 Alumina plate was chosen as the specimen. Test standards were kindly provided by both Philips and Siemens with some participants using their own standard. Participants were required to run a series of 22 scans with prescribed settings of tube current, voltage, step size, slit size and count times. One scan covered the range  $1^\circ$  to  $160^\circ$  with the others just covering a single reflection. The measurement details supplied to complete the Round Robin are given in Appendix 1. The original US Quattro Pro spreadsheet was converted into Excel 5 for UK data evaluation.

### Data Extraction

The required data was manually extracted from ASCII files submitted by participants, using Excel 5, as shown below. For tests 1, 2, 3, 5 and 6 average backgrounds from  $51.0^\circ$  to  $51.18^\circ$  and  $54.0^\circ$  to  $54.18^\circ$  were taken. The position and intensity values of the (024) alpha 1 peak and the intensities of alpha 2 peak and valley between the peaks were extracted. For test 4 the five most intense points over the alpha 1 peak and the peak position were taken. Test 7 was performed on the weak peak at  $41.7^\circ$  with a single background ( $41$  to  $41.18^\circ$ ) and the same peak information as test 1. In test 8 peak positions and intensities were taken for most peaks along with background intensities at specific angles.

## Data extraction.



## Participants.

Data for Thirty-three instruments were submitted for inclusion. A list of participants is given in Appendix 2, some included more than one instrument. Each participant was given a unique number (starting at 20) to ensure their anonymity and yet still allow them to identify their data in this report. Both Philips (19) and Siemens (13) instruments are represented with one Picker 2 circle. A range of instrument parameters was used across the instruments with both Copper and Cobalt tubes, primary and secondary beam monochromators and both fixed and variable slits. A detailed list of the usual Instrument parameters used on each instrument is given in Appendix 3.

Table 1 below summarises the instrument settings used. It shows the variety of slits used, especially for Siemens instruments, but overall it is not difficult to specify a set of generic instrument parameters as shown in bold characters. This provides a useful starting point for instrument settings for new diffraction users.

**Table 1**

<b>Instrument Statistics</b>				
		Philips	Siemens	
Model	PW1050	10		
	PW18XX	5	9	D500
	X-Pert	4	4	D5000
Tube	Co LFF	2	2	
	Cu Normal	2	2	
	Cu Broad	2		
	Cu LFF	<b>13</b>	<b>8</b>	
	Cu AEG		1	
Radius	200	1	1	
	250	2	2	
	410		2	
	Standard	<b>16</b>	<b>8</b>	
Sec Monochrom	None	2	3	
	Kevex ED		1	
	Graphite	<b>17</b>	<b>9</b>	
Geometry	Theta/Theta	1	1	
	T/2T	<b>18</b>	<b>12</b>	
Divergence Slit	Variable	6	2	
	0.3		3	
	0.5	2		
	1	<b>11</b>	<b>7</b>	
	2		1	
Receiving Slit	0.05		4	
	0.1	6		
	0.15		3	
	0.2	<b>12</b>	3	
	1	1	2	
	2		1	
Detector	Proportional	<b>19</b>		
	Scintillation		<b>12</b>	
	Kevex ED		1	

**Results.**

The exercise was divided into eight tests with each test focusing on a particular instrument parameter. The spreadsheet calculations reduced the data to a suitable figure to allow easy comparison between instruments. The calculations performed within the spreadsheet to arrive at these numbers are detailed in Appendix 4. A summary of all the calculated test results for all instruments is given in Appendix 5.

### Test 1 mA variation.

This test was used to give the comparative counts per second per mA value. The countrates at all three mA setting have been plotted against intensity for all instruments and are shown for each instrument in Appendix 6. The plots have been individually scaled for intensity so beware of direct comparison. A straight line between zero and the highest value has been added. The data should fall in straight lines if the detector is not being saturated and the dead time circuitry is functioning correctly. Non linearity with primary beam monochromators is an indication of misalignment. The plots for users 23, 38, 44, 48, and 50 are not linear and indicate problems. The test 1 data was also used to estimate deadtime. The counting statistics were not ideal to give an accurate determination. It is recommended that users with a Not Detm comment in the summary determine the deadtime more accurately for their system. Also in test 1 a Peak to background ratio (pk/Bg) was measured, a similar determination was made in test 2. For most applications a high pk/Bg is beneficial.

### Test 2 increasing kV

The increase in background with increased kV was determined. Lower values for the calculated ratio are preferred. The results are comparable with the US data. The high UK Philips average is the result of exceptionally high values obtained by users 22 and 27. An exponent value (Exp) is also calculated from the peak intensities. It should be just greater than one for copper radiation.

### Test 3 Receiving slit variation.

Participants were asked to scan using three different receiving slit sizes. A wide range of slits was used in this test and in test 5 where the divergence slit was varied. A summary of the slits employed for both tests along with the test results is shown in table 2. The intensity ratio values quoted in the summary sheet (med/sm, lg/sm) should approximate the ratio of the slit sizes used. Low values indicate that over illumination of the slit by the beam has occurred, and that the optimum slit size has been exceeded.

### Test 5 Divergence slits.

Participants were asked to use the kV and mA settings of test 3 and their usual receiving slit. Three scans were performed with 1/4, 1/2, and 1° divergence slits. The higher the reported Figure of Merit (FOM) the better the instrument performance. The UK data is comparable with USA data with Philips instruments giving larger FOM's.

Table 2

User	Test 3 Receiving Slit size			Test 3 Results		Test 5 Divergence Slit size			Test 5 Results	
	Small	Medium	Large	med/sm	lg/sm	Small	Medium	Large	Small	Medium
21	0.1	0.2	1.57	1.60	2.51	0.17	0.50	1.00	0.88	0.14
22	0.1	0.2	0	2.04	0.00	0.25	0.50	1.00	0.73	0.51
23	0.1	0.3	3	2.81	13.49	0.30	1.00	0.00	0.87	0.14
24	0.1	0.2	4	1.90	2.78	0.25	0.50	1.00	0.99	0.17
25	0.05	0.1	1.8	2.45	5.62	0.25	0.50	1.00	0.95	0.18
26	0.1	0.2	1.8	1.73	2.36	0.17	0.50	1.00	0.84	0.16
27	0.1	0.2	0.3	1.70	2.23	0.25	0.50	1.00	0.80	0.38
28	0.05	0.1	3	1.85	4.84	0.25	0.50	1.00	1.10	0.43
29	1	3	0	0.99	0.00				0.74	0.56
30	0.018	0.05	0.6	3.12	5.04	0.30	0.60	1.00	0.88	0.18
31				1.69	2.44	0.25	0.50	1.00	0.86	0.15
32	0.018	0.05	0.6	3.36	5.48	0.25	0.50	1.00	1.01	0.23
33	0.05	0.1	6	1.89	5.19	0.25	0.50	1.00	0.77	0.27
34	0.1	0.2	3	2.27	3.46	0.10	0.30	1.00	0.88	0.16
35	0.1	0.2	2	1.97	5.82	0.60	1.00	2.00	0.29	0.42
36	0.1	0.3		2.03	No Data	0.25	0.50		1.13	0.48
37	0.018	0.05	0.6	2.74	5.39	0.10	0.30	1.00	0.68	0.49
38	0.1	0.20	0.30	1.75	2.22	0.25	0.50	1.00	0.82	0.32
39										
40	0.1	0.3		2.00	No Data	0.25	0.50		1.04	0.37
41	0.018	0.05	0.6	2.54	3.88	0.10	0.30	1.00	1.10	0.43
42	0.05	0.1	0.2	1.69	3.03	0.25	0.50	1.00	0.86	0.17
43	0.1	0.2	0.3	1.69	2.09	0.25	0.50	1.00	0.67	0.49
44	0.018	0.5	6	8.39	11.68	0.25	0.50	1.00	0.93	0.21
45				No Data	No Data				0.90	0.18
46	0.1	0.2	1	1.57	2.03	0.25	0.50	1.00	0.88	0.19
47	0.1	0.2	0.3	1.86	2.46	0.25	0.50	1.00	0.90	0.15
48	0.018	0.05	0.6	2.93	4.23	0.25	0.50	1.00	0.89	0.16
49	0.05	0.10	3.00	1.60	3.84	0.25	0.50	1.00	0.64	0.48
50	0.01	0.05	0.2	5.39	9.37	0.04	1.00	0.00		
51										
52	0.1	0.2	1	1.85	2.65	0.25	0.50	1.00	0.65	0.45
53	0.05	0.1	2	1.83	4.54	0.60	1.00	2.00	0.81	0.14
54	0.05	0.1	0.6	1.93	5.43	0.25	0.50	1.00	1.03	0.17
55	0.1	0.2	0.3	2.39	2.44	0.50	1.00	0.00	0.65	0.45

#### Test 4 Step size

Participants were asked to use the kV and mA settings of test 3 and their usual receiving slit and perform scans with 0.01, 0.02, 0.04, 0.06 and 0.08 step sizes. This test identifies the maximum step size which can be used without adversely affecting the calculated peak position. The figures quoted are the peak shifts relative to the smallest step size. It shows that 0.02 steps are safe for all instruments and 0.04 on some. Above 0.04, shifts of 0.03° to 0.04° are obtained on the 52.5° peak. The results for the UK and USA show the same trend, with little difference between manufacturers.

#### Test 6 No Secondary Soller Slit

Most participants were reluctant to remove the Soller slits. Only nine data sets are available. The data shows the significant peak intensity gain obtained (Wo/W) with only a small loss of resolution (p/v). The background (Bg) also increases in most cases with the same ratio as the peak. If low intensity is a problem then removal of the Secondary Soller slit is in most cases recommended. The UK and USA data is comparable and shows that the increase in intensity is bigger on Siemens' instruments.

### Test 7 Count time on weak peak

The effect of increasing count time per step on the weak peak at 41.7° was measured. One second and nine second times were compared for net peak intensity. The data shows that for qualitative analysis a one second counting time is adequate for detecting weak peaks. The USA and UK results are in good agreement.

### Test 8 Long scan

A scan from 1° to 160° (or as wide a range as possible) using 45kV and 35mA, 0.02 steps, 2 seconds per step and the user's usual slits was requested. The data extracted were peak positions, peak intensities and background intensities. Intensity data and peak position data from the NIST instrument in the USA were used as the standard against which all other data sets were compared. If necessary the intensity data was converted from variable to fixed slit values so that direct comparison of the intensities of all instruments could be made. Each peak intensity was scaled against the strongest peak and given a percentage relative intensity ratio. This relative intensity ratio was divided by that obtained from the NIST data set to give a ratio that should ideally be unity. This comparative data for all the peaks in the scan was averaged (Av) and the standard deviation (sd) calculated and these figures are reported in the Summary. Average intensities are lower in the UK and standard deviations are a factor of 3 worse in the UK.

Appendix 7 shows a plot of intensity ratio for all instruments. The scale on the y-axis is 0.3 to 1.9 plotted against angle on the x-axis. It is clear from this plot that many instruments show the same pattern of deviation from unity. This points to an error in the NIST intensities possibly introduced in the conversion of the spreadsheet to Excel. The peak intensity data was recalculated using the certified peak intensity data provided on the NIST SRM1976 data sheet. This removed the pattern observed in the previous data and put most instruments within the 1±0.0785 specified on the standard certificate for an instrument in control. The modified intensity ratio plots are shown in Appendix 8 and the individual intensity values are given in Table 3 below. The new intensity ratios also improve the averages and standard deviations and bring them into line with the US values as shown.

		Average	Std Dev
Philips	UK	0.86	0.31
	UK Mod	1.02	0.09
	US	0.94	0.11
Siemens	UK	0.84	0.27
	UK Mod	1.05	0.11
	US	1.05	0.14

It is now clear to see that some instruments fitted with variable divergence slits show intensity problems. The data for Cobalt tubes all show the same trend of increasing intensity deviation with angle. This is being investigated by Ron Jenkins and will be discussed at the Leeds workshop.

**Table 3 Modified Intensity Ratios**

Angle	25.58	35.16	43.36	52.56	57.5	68.2	76.88	95.24	101.08	116.56	127.68	136.06	145.2	Average	Std Dev
21	1.060	1.000	1.034	1.077	1.039	1.022	1.035	1.099	1.128	0.981	1.118	1.058	1.081	1.056	0.044
22Co	1.019	1.000	0.954	0.949	0.877	0.910	0.969	1.229	1.279					1.021	0.139
23	0.969	1.000	1.108	1.142	1.032	0.977	1.030	0.987	1.058	0.851	0.992	1.045	1.082	1.021	0.073
24	1.093	1.000	1.007	0.999	0.955	1.005	0.978	1.145	1.067	1.197	1.145	1.231	1.211	1.079	0.097
25	0.875	1.000	0.905	0.872	0.886	0.919	0.980	1.027	1.046	0.958	0.995	0.982	1.053	0.961	0.064
26	1.002	1.000	0.966	0.944	1.036	0.964	0.988	1.016	1.009	1.013	0.988	1.046		0.998	0.030
27	0.954	1.000	1.009	1.059	1.024	1.065	1.093	1.142	1.131	1.172	1.122			1.070	0.068
28v	1.061	1.107	1.105	1.034	1.000	0.962	1.019	0.924	0.944	0.974	0.863	0.819	0.873	0.976	0.091
29	1.153	1.000	0.949	1.019	0.941	0.998	0.616	1.499	1.755					1.103	0.336
30	0.955	1.000	1.018	1.039	1.058	1.057	1.123	1.081	1.147	1.149	1.082	1.029	1.002	1.057	0.059
31	1.040	1.000	1.030	1.045	1.049	0.912	1.009	1.070	1.019	1.096	1.083	1.036	1.102	1.038	0.049
32	0.977	1.000	0.956	0.979	0.978	1.028	1.117	1.154	1.166	1.213	1.286	1.300	1.478	1.126	0.161
33v	1.061	1.000	0.966	0.959	1.020	1.003	0.965	0.933	0.934	1.022	0.949	0.948		0.980	0.040
34	1.006	1.000	1.016	1.020	1.041	1.129	1.065	1.133	1.070	1.064	1.112	0.986	1.151	1.061	0.056
35	1.035	1.000	1.011	1.074	1.052									1.034	0.030
36v	1.144	1.025	1.009	1.050	1.000	0.889	1.046	0.990	0.995	0.977	0.943	0.944	0.984	1.000	0.062
37Co	1.075	1.000	0.924	0.966	0.991	0.936	1.098	1.278	1.347					1.068	0.151
38v	1.365	1.285	1.178	1.101	1.000	1.021	0.854	0.717	0.719	0.608	0.566	0.531		0.912	0.287
40v	1.139	1.078	1.004	1.005	1.000	0.988	0.939	0.904	0.922	0.726	0.866	0.797	0.872	0.942	0.112
41	0.995	1.000	0.972	1.015	0.989	0.996	1.069	1.157	1.195	1.276	1.325	1.260	1.427	1.129	0.154
42	0.873	1.000	0.924	0.888	0.994	0.982	1.039	1.007	0.940	0.979	0.949	0.992	1.030	0.969	0.051
43vCo	0.852	1.000	1.033	1.060	1.026	1.056	1.119	1.278	1.476					1.100	0.180
44	0.958	1.000	1.014	0.803	1.023	0.929	1.099	1.285	1.126	1.181	1.198	1.206	1.224	1.081	0.140
45	1.006	1.000	1.010	1.067	1.044	1.114	1.089	1.053	1.092	1.109	1.034	1.057	1.106	1.060	0.040
46	1.006	1.000	1.010	1.067	1.044	1.114	1.089							1.047	0.045
47	1.039	1.000	0.990	0.996	0.969	0.873	0.976	1.056	0.943	1.003	1.006	1.030	1.119	1.000	0.058
48	0.995	1.000	1.053	0.984	1.027	0.938	1.086	1.122	1.126	1.091	1.121	1.095	1.013	1.050	0.061
49v	1.320	1.280	1.070	1.070	1.000	1.010	0.900	0.860	0.770					1.031	0.182
50	0.951	1.000	1.046	1.046	1.044	1.155	0.952	0.976	1.062	0.842	1.072	1.148	1.013	1.024	0.084
52	0.965	1.000	0.954	0.965	0.971	1.002	1.044	1.151	1.072					1.014	0.065
53	1.050	1.000	0.963	0.963	1.026	0.872	0.884	0.943	0.951	0.910	0.952	0.966	0.975	0.958	0.051
54v	1.106	1.069	0.999	0.975	1.000	1.085	0.998	0.927	0.893	0.978	0.856			0.990	0.078
55	1.042	1.000	0.977	0.960	0.980	1.056	0.992	1.032	1.076					1.013	0.040

Average Siemens	1.051	0.107
Average Philips	1.015	0.092

Appendix 8 shows the background data for all instruments. The scales are; intensity ratio to the NIST instrument for the y-axis and angle for the x-axis. Even a range of 0 to 3.5 in the ratio is not wide enough to fit all the data from some instruments. The data from 22,35 and 47 show the high background values at low angle expected from the use of a beta filter monochromator. The data from 29 shows a peak at 10° with a 7.2 ratio. The data from 21,27 and 52 show high values of background intensity ratio at 5° with 19, 27 and 29 ratios respectively, possibly due to sample holder scatter.

Each participant was provided with a plot showing peak displacement relative to the NIST data. The plots show that many instruments have displacement problems. These plots are shown in Appendix 10. They all are on the same scale of displacement -0.12° to +0.12° on the y-axis and angle on the x-axis. The user number is given at the top of each plot, the addition of “Co” signifies data from a Cobalt tube and “v” that a variable divergence slit was employed. The data shows that many instruments have angular displacement problems. Some have a simple displacement over the whole angular range e.g. user 40 which is easily corrected. Others have displacement that varies with angle e.g. user 29 which is symptomatic of a one to two theta alignment error.

### **Instrument Benchmark**


From the US data Jenkins and Schreiner derived a Figure of Merit (FOM) to test instrument sensitivity - it is fully explained in a book by Jenkins and Snyder(ref 1). I have taken Test 1c

UK data to calculate the FOM. this gives participants a benchmark to test their instrument sensitivity against other similar instruments

The data for 10 and 46 are from the same instrument with a change in X-ray tube, all other parameters remained the same. The 46 data set was obtained with a new X-ray tube, of the same specification ( Cu LFF), as an old one used for the 10 data.


User ID	FOM
37	4.3
23	4.4
50	4.8
53	6.1
10	7.5
33	7.8
34	8.5
44	8.6
22	8.9
41	9.0
48	9.1
28	9.5
43	9.5
45	10.8
52	12.0
49	12.5
36	12.8
40	12.9
47	13.1
54	13.2
29	13.5
42	13.9
30	14.2
55	14.4
21	15.0
25	15.0
31	15.2
24	15.7
26	15.9
38	18.1
35	19.6
32	21.2
27	22.5
46	26.8

Figure of Merit



## Calculating the FOM

### Instrument Parameter Figure of Merit



$$FOM = \frac{MW}{\sqrt{W(M + 4B)}}$$

- M = Net Intensity cps
- B = Average Background
- W = FWHM

The higher the value of the FOM the better the sensitivity. The FOM's derived form a useful basis for testing the performance of any instrument. The value calculated for an individual instrument can be used to validate that any changes in alignment or measurement parameters have improved performance.

## Conclusions

There is a wide spread of performance, many instruments have angular displacement problems which are easily overcome. The errors found demonstrate the need to frequently check goniometer angles and calibrate at more than just the single angular position which is common practice amongst the diffraction community. Count rates are significantly lower than



in the US. This could be due to X-ray tube degradation. Perhaps routine instrument checks should be more widely used to ensure that equipment performs to manufacturer's specification and that X-ray tube degradation is monitored. Routine checks should monitor overall or angular dependent loss in intensity or shifts in 2 theta values and changes in peak shape. The NIST SRM 1976 is a useful standard to use for such tests especially where they are performed as part of an accreditation procedure. Test 8 could be used as the basis for the test but modified to use your normal power setting. These data could then be used to generate the displacement, background and intensity plots and the 52.5° peak could be used to give the instrument sensitivity FOM. All instrument manufacturers should supply an instrument routine check programme and certified standard with their systems.

## Reference

1. *Introduction to X-Ray Powder Diffractometry - Jenkins & Snyder - John Wiley & Sons Inc.*  
p224-226 ISBN 0-471-51339-3

## APPENDIX 1

### Measurement instructions UK Sensitivity Round Robin

**USER NUMBER.....** a unique identifier supplied by the co-ordinator

**Serial Number of SRM 1976 used .....**(on front of standard box)

All measurements to be made using NIST SRM1976 Alumina plate

All angles quoted are for a Copper tube and will require conversion for other anode materials.

Where possible the following file name convention should be used:-

**\*\*uk??** ( Where \*\* is User identity number and ?? is test number e.g. 1a)

**For tests 1 to 5 the following conditions apply:-**

Scan 50-55° , Step size 0.02 ° (except test 4), count time 4 seconds/step

**MAKE SURE YOU DO NOT EXCEED YOUR TUBE LOADING.**

**Test 1** Calculate X from  $X = \text{Max. tube loading in watts} \div 180$

e.g. for 2 kW tube  $X = 2000 \div 180 = 11.11$

X is rounded to a mA interval suitable for your

generator (e.g. for 11.11 use 10 or 11 mA)

Use your usual slits and 45 KV for this test

**Divergence slit used.....\*Receiving slit used.....Detector type.....eg Scint**

**\*Mark VAR for variable divergence slit**

1a Scan at 1X mA (eg 1.8 kW tube use 10 mA) **mA used.....File.....**

1b Scan at 2X mA (eg 1.8 kW tube use 20 mA) **mA used.....File.....**

1c Scan at 4X mA (eg 1.8 kW tube use 40 mA) **mA used.....File.....**

**Test 2** Use your usual slits and 20 mA

2a Scan at 30 KV **File name.....**

2b Scan at 40 KV **File name.....**

2c Scan at 45 KV **File name.....**

2d Scan at 50 KV **File name.....**

**Test 3** Use your usual divergence slit, 45 KV and 35 mA (or as close as your generator will allow). **Values used KV..... mA.....**

3a Scan using smallest available receiving slit **Value used.....File.....**

3b Scan using next smallest avail receiving slit **Value used.....File.....**

3c Scan using largest available receiving slit **Value used.....File.....**

**Test 4** Use your usual slits and Test 3 power

4a Use 0.01° step size **File name.....**

4b Use 0.02° step size **File name.....**

4c Use 0.04° step size **File name.....**

4d Use 0.06° step size **File name.....**

4e Use 0.08° step size **File name.....**

- Test 5**                      **For fixed divergence slit systems only**  
 Use usual receiving slit and Test 3 power
- 5a                      Scan with 1/4° divergence slit                      *File name.....*  
 5b                      Scan with 1/2° divergence slit                      *File name.....*  
 5c                      Scan with 1° divergence slit                      *File name.....*
- Test 6**                      Remove secondary soller collimator (if possible) and repeat  
 test 4b                      *File name.....*
- Test 7**                      Obtain profile of weak peak at 41.7° by scanning 41° to 42.5°  
 step size 0.02° , your usual slits and Test 3 power
- 7a                      Scan at 1 second/step                      *File name.....*  
 7b                      Scan at 9 second/step                      *File name.....*
- Test 8**                      Perform as long a scan as possible starting at 1° 2 Θ and  
 scanning to say 160° 2 Θ. Use Test 3 power, 0.02° step size,  
 2 seconds / step and your usual slits.                      *File name.....*

**Reporting results - For all tests:-**

- 1                      Plot each scan with the largest peak at about full scale For test 8 split into 20° segments and scale individually
- 2                      Computer printout of raw count data vs. 2Θ (Disable any automatic angular correction made by your system) or raw data onto IBM PC compatible floppy

**APPENDIX 2**

**UK Round Robin Participants**                      *Italics Siemens / Normal- Philips + Picker 2 circle*

*D Beveridge - Ilford Ltd*  
 Mr Bruce Fox - Raychem Ltd  
*Ms Michele Leigh- De Beers Industrial Diamond Division*  
 Mr NL Andrew - Philips Research Laboratories  
*Mr R R Giles - Alcan Chemicals Ltd*  
 Dr Andries Bosland /Mr Jes Brown - Alcan International Ltd  
 Dr Craig Adam - Staffordshire University  
*Mr F L Cullen - AEA Technology*  
 Mr Eric Kelly - 44 Penymynydd Rd  
 Dr J M Parker - Dept of Engineering Materials Sheffield University  
 Dr Mary Vickers - Cambridge University  
*Dr Mary Vickers - Cambridge University*  
 Dr Mary Vickers - Cambridge University  
*Dr Tim Hyde - Johnson Matthey Technology Centre*  
*Dr Robin Payne - Physical Sciences Group Zeneca*  
*Mr Mike Edmondson - Research Dept Zeneca*  
*Mr Mike Edmondson - Research Dept Zeneca*  
*Mr Ian Tucker Sharon Evans - Unilever Research*  
*Dr J A Jutson- BICC Cables Ltd*  
 Mr T J Carter - AWE Aldermaston  
 Mrs J Blake- ICI C&P Ltd Runcorn  
 Dr David Dyson - British Steel Technical Rotherham  
 Mr David L Rush - Hepworth Refractories Ltd  
 Mr David J Taylor - Pilkington  
*Dr S V Norval- ICI Wilton*  
*Ms Cerys Oliver - Cookson Technology Centre*  
 Miss Sue Staniforth - GalaxoWellcome  
 Mr Colin Small - Rolls-Royce PLC  
 Mr P G Lake - Glaxo Wellcome  
 Dr Christian Lehmann - Chemistry University of Durham  
 Dr R Alan Howie - University of Aberdeen  
 Mr A Calder - Dept of Geology St Andrews University  
 Dr J Ian Langford - School of Physics Birmingham

**APPENDIX 3**

See Excel spreadsheet rrobpnd.xls after appendix 4

## **APPENDIX 4**

### **Spreadsheet Calculations**

#### **Test1**

**c/s/mA** - This is the background subtracted, dead time corrected, counts per second per milliamp.

**Dead Time** - This is the counter dead time calculated from the peak intensities at different mA settings. Note that the counting statistics are not ideal for the calculation. If your data looks odd you should recalculate with better statistics.

**Pk/Bk** - This is calculated from test 1b as  $(I-B)/B$   
where I = peak Intensity (cps) B = Ave Background (cps) . The background is averaged from 51.0 to 51.18 and 54.0 to 54.18 degrees.

#### **Test2**

**Bk inc with kV** - This is  $((A-B)/B)/(C-D)$   
where A is Ave background/kV from 45 kV data  
where B is Ave background/kV from 30 kV data  
Where C is kV -8.86 for 45kV data  
Where D is kV -8.86 for 30kV data

**Pk/Bk (40kV)** - This  $((A-B)/40)/(B/40)$   
where A is peak intensity from 40 kV expt  
where B is average background intensity from 40 kV expt.

**Exponent Value** - This is  $(I1-I2)/(k1-k2)$   
where I1 is log peak intensity from 45kV  
where I2 is log peak intensity from 40kV  
where k1 is log (45-8.86)  
where k2 is log (40-8.86)

#### **Test3**

**Int Ratio x/y** this is the Peak-Bk intensity of slit x divided by the Peak-Bk of slit y

#### **Test4**

**Pk Shift 0.01/0.02** - The alpha 1 peak intensity = T with the adjacent intensities on the low angle side = R and S and on the high angle side U and V all in counts per second.

Pk Shift = A-B

where A =  $C-(D/E)*F$  for 0.02 step measurements

C is peak position as 2 theta angle

D is  $-0.2R-0.1S+0.1U+0.2V$  (1st derivative)

E is  $0.285714R-0.142857S-0.285714T-0.14857U+0.285714V$  (2nd derivative)

F is step size

B =  $C-(D/E)*F$  for 0.01 step measurements

The other Pk Shift data is calculated as above but substituting the 0.04 and 0.08 for the 0.02 step data.

#### **Test5**

**FOM Small** -

$$A1 * \sqrt{(1/(A2/V)/(A1+4*B))}$$

where A1 is alpha1 net peak intensity for small divergence slit

where A2 is alpha2 intensity

where V is valley intensity

where B is the average background intensity

**FOM Mid** - as above but using middle divergence slit data

#### **Test6**

**Intensity gain Wo/** - c/s/mA without collimator divided by c/s/mA with collimator

**Change Peak/Valley** - (A/V)/(B/W)

where A is alpha2 intensity without collimator

where V is valley intensity without collimator

where B is alpha2 intensity with collimator

where W is valley intensity with collimator

**Change Bk** - average background without divided by average background with sec collimator

#### **Test7**

**Improvement factor** - (A/B)/ $\sqrt{t}$

where t is 9 second count time

A is (I\* $\sqrt{t}/\sqrt{C}$ ) for 9 second data

where I is net peak intensity

where t is count time

where C is average background

B is (I\* $\sqrt{t}/\sqrt{C}$ ) for 1 second data

#### **Test8**

Peak Intensities were corrected :-  $I_{cor} = I * \sin(\text{rad}(A/2))$  and scaled to give  $I/I_{max} * 100$

where A is the peak position as 2 theta angle and I is peak intensity.

Your relative intensities were then ratioed to the intensities obtained on the NIST instrument .

The **Average intensity** reported is the mean of all your individual test 8 relative intensity ratios and the

**Standard deviation** is the 1 sigma value of this intensity data.

#### **Appendix 3 & 5 to 10**

These pages are from the Excel spreadsheet rrobapnd.xls and follow on the next page of this PDF version.

ApX 3 Instrm

ID	Manufact	Model	Anode	Filament	TakeOff	Radius	PB mono	SB moNone	Geometry	Div	Rec	Counter
21	hybrid	PW1050	Cu	Normal	6	Standard		graphite	T/2T	1	0.2	prop
22	Philips	PW1820	Cobalt	LFF	6	Standard		None	T/2T	1	0.1	prop
23	Siemens	D500	Cu	Normal	6	250	Quartz	None	T/2T	1	1	Scint
24	Philips	PW1050	Cu	Broad	6	Standard		graphite	T/2T	1	1	prop
25	Philips	X'pt PW3020	Cu	LFF	6	Standard		Graphite	T/2T	0.5	0.2	prop
26	Philips	PW1820	Cu	LFF	6	250		graphite	T/2T	1	0.2	prop
27	Philips	PW1050	Cu	LFF	6	Standard		graphite	T/2T	1	0.2	prop
28	Philips	X-Pert MPD	Cu	LFF	6	200		graphite	Theta/Theta	var	0.2	prop
29	Siemens	D500	Cobalt	FF	6	Standard		None	T/2T	1	0.05	Scint
30	Siemens	D500	Cu	LFF	6	Standard		graphite	T/2T	1	0.05	Scint
31	Philips	X-pert	Cu	LFF	6	Standard		graphite	T/2T	1	0.2	prop
32	Siemens	D500	Cu	Normal	6	250		graphite	T/2T	1	0.15	Scint
33	Siemens	D5000	Cu	LFF	6	200.5		graphite	Theta/Theta	var	0.2	Scint
34	Siemens	D500	Cu	LFF	6	Standard		graphite	T/2T	0.3	1	Scint
35	Siemens	D5000	Cu AEG	fk600412	6	414		None	T/2T	2	2	Scint
36	Philips	PW1800	Cu	LFF	6	Standard		graphite	T/2T	var	0.1	prop
37	Siemens	D500	Cobalt	LFF	6	Standard		graphite	T/2T	1	0.05	Scint
38	Philips	PW1050	Cu	LFF	6	Standard		graphite	T/2T	var	0.2	prop
39	Philips	PW1050	Cu	Broad	6	Standard		graphite	T/2T	No DATA Supplied		
40	Philips	PW1800	Cu	LFF	6	Standard		graphite	T/2T	var	0.1	prop
41	Siemens	D500	Cu	FF	6	Standard		graphite	T/2T	0.3	0.15	Scint
42	Philips	PW1050	Cu	LFF	6	Standard		graphite	T/2T	1	0.1	prop
43	Philips	PW1820	Cobalt	LFF	6	Standard		graphite	T/2T	var	0.1	prop
44	Siemens	D500	Cu	LFF	6	Standard		graphite	T/2T	0.3	0.15	Scint
45	Philips	PW1050	Cu	LFF	6	Standard		graphite	T/2T	0.5	0.2	prop
46	Phillips	1050mod	Cu	LFF	6	Standard		graphite	T/2T	1	0.2	prop
47	Philips	PW1050/80	Cu	Normal	6	Standard		None	T/2T	1	0.1	prop
48	Siemens	D500	Cu	LFF	6	410		graphite	T/2T	1	0.05	Scint
49	Philips	PW3040	Cu	LFF	6	250		graphite	T/2T	var	0.2	prop
50	Picker	2-Circle	Cu	LFF	3	145.5	Ge 111	None	T/2T	2	0.05	scint
52	Philips	PW1050	Cu	LFF	6	Standard		graphite	T/2T	1	0.2	prop
53	Siemens	D5000	Cu	LFF	6	Standard		graphite	T/2T	1	0.2	Scint
54	Siemens	D5000	Cu	LFF	6	Standard		ED -Det	T/2T	var	0.2	ED KeveX
55	Philips	PW1050/81	Cu	Broad	6	Standard		graphite	T/2T	1	0.2	prop

UK Intensity Round Robin  
Data Summary

Appendix 5 ID	Type	Test #1 c/s	#1 DT	#1 PkBg	#2 Bg/kV	#2 Pk?bg	#2 Exp	#3 med/sm	#3 Ig/sm	#4 1/2	#4 1/4	#4 1/8	#5 med/sm	#5 Ig/sm	#6 Wo/W	#6 P/V	#6 P/B	#7 IF	#8 Av	#8 SD
21 Philips		58	22	58	1.32	64	0.83	1.60	2.51	0.0019	0.0012	0.0386	12.50	21.05	No Data	No Data	No Data	1.00	0.88	0.14
22 Philips Co		26	15	27	7.28	37	1.03	2.04	0.00	-0.0167	-0.0184	-0.0212	No Data	No Data	No Data	No Data	No Data	1.37	1.02	0.23
24 Philips		42	Not Detm	47	1.84	47	0.97	1.90	2.78	0.0071	0.0296	0.0407	18.32	26.62	No Data	No Data	No Data	1.06	0.99	0.17
25 Philips		62	11	67	1.53	65	1.01	2.45	5.62	0.0011	0.0044	0.0338	21.10	34.62	No Data	No Data	No Data	1.01	0.95	0.18
26 Philips		82	Not Detm	69	1.19	68	1.03	1.73	2.36	-0.0017	-0.0068	0.0356	14.32	28.58	1.78	0.75	2.00	0.91	0.84	0.16
27 Philips/hybrid		109	11	61	1.58	59	1.20	1.70	2.23	0.0075	-0.0053	-0.0393	No data	32.49	No Data	No Data	No Data	1.06	0.93	0.19
28 Philips		34	Not Detm	82	3.07	73	1.42	1.85	4.84	0.0077	0.0186	0.0428	No data	No data	1.95	0.88	2.24	1.02	0.85	0.17
31 Philips		63	5	69	1.71	67	1.17	1.69	2.44	0.0011	0.0270	0.0496	13.70	25.94	1.80	0.83	1.96	0.95	0.86	0.15
36 Philips		53	Not Detm	83	1.07	87	1.19	2.03	No Data	-0.0008	0.0009	0.0342	No Data	No Data	No Data	No Data	No Data	0.96	0.86	0.17
38 Philips		57	Not Detm	49	1.24	62	0.86	1.75	2.22	-0.0092	-0.0158	0.0149	No Data	No Data	No Data	No Data	No Data	0.93	0.78	0.24
40 Philips		50	7	84	3.31	80	1.09	2.00	No Data	-0.0023	-0.0020	0.0523	No Data	No Data	No Data	No Data	No Data	1.01	0.81	0.15
42 Philips		58	10	67	1.17	66	1.35	1.69	3.03	-0.0032	0.0037	0.0316	15.97	23.42	No Data	No data	No Data	1.14	0.86	0.17
43 Philips Co		30	6	107	0.41	110	0.61	1.69	2.09	-0.0028	-0.0086	0.0337	No Data	No Data	No Data	No Data	No Data	1.04	0.80	0.25
45 Philips		No Data	No Data	No Data	No Data	No Data	No Data	No Data	No Data	No Data	No Data	No Data	No Data	No Data	No Data	No Data	No Data	No Data	0.90	0.18
46 Philips		181	5	69	1.26	68	0.90	1.57	2.03	-0.0025	-0.0103	0.0155	31.69	45.41	No Data	No Data	No Data	0.98	0.88	0.19
47 Philips		57	13	32	5.39	38	1.14	1.86	2.46	0.0061	0.0148	0.0370	14.87	20.44	No Data	No Data	No Data	1.06	0.90	0.15
49 Philips		43	Not Detm	87	1.84	94	0.88	1.60	3.84	-0.0003	0.0011	0.0278	No Data	No Data	No Data	No Data	No Data	1.08	0.89	0.28
52 Philips		45	<10	47	No Data	48	No Data	1.85	2.65	0.0066	-0.0122	0.0002	15.29	21.67	No Data	No Data	No Data	0.98	0.91	0.19
55 Philips		34	8	62	1.50	64	1.26	2.39	2.44	-0.0100	-0.0074	0.0224	15.69	20.78	1.92	0.88	2.19	1.02	0.91	0.16
30 Siemens		52	14	69	1.39	73	0.83	3.12	5.04	0.0020	-0.0088	0.0327	22.12	28.67	2.90	0.65	4.08	1.07	0.88	0.18
48 Siemens		40	Not Detm	75	2.69	76	1.82	2.93	4.23	0.0005	-0.0061	0.0363	13.90	24.71	No Data	No Data	No Data	0.99	0.89	0.16
54 Siemens		57	17	81	1.45	90	0.73	1.93	5.43	0.0092	0.0207	0.0470	9.10	17.29	2.76	0.61	3.23	1.09	0.86	0.17
35 Siemens		136	4	14	0.87	15	1.18	1.97	5.82	-0.0028	0.0313	0.0256	25.73	32.30	No soller used in all tests			No Data	0.82	0.13
33 Siemens		24	Not Detm	64	2.47	72	0.68	1.89	5.19	0.0007	0.0072	0.1721	9.23	12.63	2.89	0.76	1.55	0.92	0.84	0.16
44 Siemens		12	Not Detm	51	0.17	56	1.44	8.39	11.68	-0.0006	0.0115	0.0330	19.00	26.33	2.57	0.91	3.21	0.89	0.93	0.21
32 Siemens		80	14	35	1.05	34	0.95	3.36	5.48	0.0001	0.0106	0.0208	15.61	28.08	No Data	No Data	No Data	0.96	1.01	0.23
41 Siemens		15	Not Detm	45	2.03	48	1.03	2.54	3.88	-0.0045	-0.0031	0.0343	2.88	21.44	No Data	No Data	No Data	0.88	1.10	0.43
34 Siemens		20	Not Detm	73	1.79	74	0.69	2.27	3.46	0.0000	-0.0074	0.0343	9.57	17.48	No Data	No Data	No Data	0.98	0.88	0.16
29 Siemens Co		80	Not Detm	30	3.25	32	1.10	0.99	0.00	-0.0023	-0.0017	0.0352	No Data	No Data	No Data	No Data	No Data	1.13	1.03	0.33
37 Siemens Co		17	Not Detm	15	0.51	82	1.44	2.74	5.39	-0.0004	0.0095	0.0362	3.23	5.19	No Data	No Data	No Data	1.07	0.95	0.24
23 Siemens pbm		2	Not Detm	99	6.24	92	2.83	2.81	13.49	0.0014	0.0008	0.0030	No Alpha2	No Alpha2	No Data	No Data	No Data	0.87	0.87	0.14
53 Siemens		11	Not Detm	78	0.67	66	1.09	1.83	4.54	-0.0002	0.0008	0.0434	7.71	9.14	2.76	0.72	2.66	0.63	0.81	0.14
50 Picker pbm		4	Not Detm	75	-1.73	59	1.35	5.39	9.37	-0.0060	-0.0100	0.0120	No Alpha2	No Alpha2	4.22	No Alpha2	6.20	0.82	0.85	0.13
Average Philips UK		60	10	65	2.15	67	1.06	1.86	2.72	-0.0006	0.0008	0.0250	17.34	27.37	1.86	0.84	2.10	1.03	0.89	0.18
Average Philips USA		148	13	72	1.57	62	1.03	1.93	2.61	-0.0004	0.0023	0.0447	20.65	32.20	1.68	0.89	1.69	1.01	0.94	0.11
Average Siemens UK		42	12	56	1.89	62	1.22	2.83	5.66	0.0002	0.0050	0.0426	12.55	20.30	2.78	0.73	2.95	0.96	0.91	0.21
Average Siemens USA		56	7	72	2.00	72	1.26	1.67	2.99	0.0020	0.0049	0.0381	10.51	18.35	3.20	0.99	4.08	0.94	1.05	0.14

