


Validation and checking of crystal structures

refine_ls_number_reflns	6811
refine_ls_number_parameters	261
refine_ls_number_restraints	0
refine_ls_R_factor_all	0.0455
refine_ls_R_factor_gt	0.0284
refine_ls_wR_factor_ref	0.0560
refine_ls_wR_factor_gt	0.0537
refine_ls_goodness_of_fit_ref	0.923
refine_ls_restrained_S_all	0.923
refine_ls_shift/su_max	0.002
refine_ls_shift/su_mean	0.000



1



OUTLINE

- Overview validation and checking
- Validation for Acta C, etc.
- Validation for other journals
- The limits of validation

3

Validation involves comparison against a set of test criteria

- Do cell volume and cell parameters match?
 - Do bonded atoms have compatible U^{eq} values?
 - Has the refinement converged?
 - Is the space group correct?
 - Are the assigned atom types correct?
- etc, etc, etc

4

Valid-ation

Correct

Appropriate

Defensible

Checking is additional to validation

- Does the structure make sense to you?
- Does the structure look right?
- Do chemically equivalent bonds agree?
- Are all CIF entries complete and correct?

6

Automated data validation with checkCIF or PLATON

- Checks for
 - CIF construction and syntax errors
 - missing information
 - parameters outside expected norms
 - conformation with convention

7

ALERT LEVELS

A Serious – attention essential

Item omitted or large deviation from norm

Alert A No crystal dimensions have been given

Alert A Ratio of Tmax/Tmin expected is > 1.30
An absorption correction is required.

Alert A Atom C58A ADP max/min Ratio 18.00

8

ALERT LEVELS

B Significant – action needed?

Item is a significant or unexpected outlier

Alert B The formula has elements in wrong order

Alert B ADDSYM detects Cc to Fdd2 transformation

Alert B Refined extinction parameter < 1.9s

Alert B Structure contains VOIDS of 130.00 Å³

9

ALERT LEVELS

C Outside expected norms – examine

May appear trivial, but do not dismiss out of hand
- an extensive list may indicate problems

Alert C Moiety formula not given

Alert C Short inter X...Y contact: O7...C1 = 2.96 Å

Alert C Low U(eq) as compared to neighbors: C1

Alert C D-H without acceptor N2-H2

C1 and N2 should be N and C, respectively

10

A/B/C indicate the seriousness of the problem

ALERT Type 1:
CIF construction/syntax error, inconsistent or missing data

ALERT Type 2:
Indicator that the structure model may be wrong/deficient

ALERT Type 3:
Indicator that the structure quality may be low

ALERT Type 4:
Cosmetic improvement, query or suggestion

Not all combinations are logical, for example 4 A

11

Sources of outlier parameters

- Unresolved feature (e.g., untreated disorder)
- Artefact due to limited data quality
- Inadequate procedures (e.g., poor corrections)
- Incorrect structure (e.g., wrong space group)
- A genuinely unusual observation!!

12

What does validation software do?

- Identifies possible problems *via* ALERTs
- Provides explanations of ALERTs
- Suggests interpretations and possible solutions

Not just for authors


- referees use it for assessment
- authors need to be aware of this
- how appropriate are IUCr criteria?

When to validate?

- software for data collection, refinement, etc
 - *should* do its own validation
- use PLATON in final stages of determination
- validate raw CIF from the refinement program
- must validate the *final* version as well
- avoids problems at submission, refereeing, etc



14

Looking at the structure



A visual examination can often be revealing: here there are some extreme ellipsoids which are also incompatible with a rigid bond model


15



P.J. Cox, RGU, Aberdeen

17

Less satisfactory



Anon

18

VALIDATION/CHECKING PROCEDURE

1. Check the CIF from refinement using PLATON
2. Augment CIF using e.g. XCIF and enCIFer
3. Re-check the CIF using PLATON or checkCIF
4. Look at ellipsoid plots from several directions
5. Check bond lengths are sensible and consistent
6. After any changes, re-check the CIF

19

Validation and IUCr Journals



Early 1990's - CIF introduced

- allows automatic creation of tables
- enables full electronic submission/processing
- increases efficiency, faster publication times
- automates many editorial tasks
- improves appearance of the journal

→ permits automated validation ←

21

Automation of syntax and data checks

- authors get instant, anonymous feedback
 - can detect and fix problems before submission
 - fewer, shorter revision cycles
 - consistent application of acceptance criteria
 - editors/referees can focus on science
-
- **RESULT: faster publication times**

22

How to get a CIF through

- Give ALL Alerts due consideration
 - appreciate validation criteria
 - criteria are based on normally expected results from routine analyses
 - Why, then, is your structure not routine?
- In any VRF...
 - avoid casual or circular responses
 - show you understand the causes of the outlier
 - explain why it is a true feature of the analysis


23

Common problems ...

- ✗ Data completeness or resolution too low
- ✗ Maltreatment of H atoms
- ✗ Structure not at convergence
- ✗ Missing or inadequate absorption correction
- ✗ Indications of a poor structure

24

All Acta C submissions in 2000




25


checkCIF in 2005

- the new home of checkCIF: <http://checkcif.iucr.org>
- service sponsored by ACS, CCDC and Elsevier
- an ORTEP plot is now included
- part of new Acta C/E submission procedures
- online upload of all material for Acta C and E papers (CIF + figures/schemes/structure factors)
- coming in the next 2-3 months

26



27



28

Validation and other Journals

29

Standards, procedures vary widely

- some journals perform extensive checks
 - some do only very basic checks
 - some do none at all
- ? so what do authors do ?

Perform your own validation

- ✓ ensure there are no serious mistakes
- ✓ ensure the quality is adequate
- ✓ submit a copy of the checking report

30

Example – an ACS procedure

- Authors submit the CIF along with the paper
- CIF must contain author names and paper title
- authors must have checked the CIF first
 - the check report may be requested
- reviewers have Web access to the CIF, along with the manuscript and any supplementary data

31

Other procedures

- Submit paper to journal
 - get a code for the paper
 - submit CIF under this code
- Submit CIF to CCDC or ICSD
 - get deposition number
 - include number in paper

32

If you have a “difficult” structure

- Identify and describe the problem
- Give details of the remedial action taken
- Describe the (successful?) outcome

33

Where and how?

1. Briefly in any experimental footnote
2. At the top of the CIF
 - use `_refine_special_details`
3. In any other Supplementary Data

34

Example of text

`_refine_special_details`

Disorder was identified in one of the tetrafluoroborate anions. All the F atoms were affected and two orientations were identified.

Similarity restraints were applied to B-F distances, and to F-B-F angles. All F atoms were refined isotropically. The occupancies of each group of four partially-occupied F atoms were refined competitively using a free variable. Each F atom was found to be disordered over two equally occupied sites, as shown by the final group occupancies of 0.506(12) and 0.494(12).

In the final model the range of B-F distances was 1.31(2)-1.42(2) Angstroms and the F-B-F angles spanned 105(2)-112(2) degrees. No difference Fourier peak in the region exceeds 0.6 e/A³.

;

35

The limits of validation

(automatic validation will not catch every problem)

36


Possible limits to validation

- × test not (yet) implemented
- × test not practical
- × error not a validation issue
- × error cannot be detected from data in CIF
- × nonsense entries in the CIF

37

Test not implemented

Example: High ADPs on isolated atom




Not detected by

- × rigid bond test
- × atom type test
- × ADP ratio test

Atom is probably O rather than Cl

38

Test not practical



C-C range is 1.49 to 1.60 Å


However, C-C single bonds are found within this range

39

Error not a validation issue

<code>_chemical_formula_sum</code>	<code>'C24 H12 Fe O6'</code>
<code>_exptl_crystal_description</code>	needle
<code>_exptl_crystal_colour</code>	colourless
<code>_exptl_crystal_size_max</code>	0.28
<code>_exptl_crystal_size_mid</code>	0.24
<code>_exptl_crystal_size_min</code>	0.03

Error not detectable from CIF data




Prior chemical information:
Complex is either Ru/Ru or Ru/Zn


Refinement as Ru/Ru gave $R_1 = 0.064$; unusual five-coordinate Ru geometry

41


Difference map with Ru/Ru model ($R_1 = 0.064$)



Difference map with Zn/Zn model ($R_1 = 0.022$)



Other examples




Ag- **CN**-Ag link,
Ag- **NC**-Ag link
or disordered?

could only be
resolved
using **DF maps**

44

Model complexes for [NiFe] hydrogenase



need good data to
distinguish Ni and Fe by
refinement

Ni ($Z = 28$) vs Fe ($Z = 26$)


$\Delta Z/Z = 2/28$

$R_1 < 4 \Delta Z?$

Use **DF maps** to confirm
and when data are poor

45

Lanthanide complexes



Ln = Er, Tm or Yb ?

$Z = 68, 69, 70$

Similar co-ordination

Similar geometry
parameters

Crystallography is not much good at distinguishing these metals.

46

Nonsense entries in the CIF *

```
_diffrn_ambient_temperature      293(2)
_diffrn_radiation_wavelength     0.69010
_diffrn_radiation_type           synchrotron
_diffrn_radiation_source         'fine-focus sealed tube'
_diffrn_radiation_monochromator  'graphite'
_diffrn_measurement_device_type  'SMART 1k on Daresbury SRS Station 9.8'
```

* see W. Clegg, *Acta Cryst.* 2003, **E59**, e2-e5

47

Unsuitable SHELX(T)L-97 defaults ?

- space group notation
- diffractometer
- $T = 293 \text{ K}?$
- absorption correction
- total data collected
- index limits
- R_{int}
- H atom treatment
- weighting scheme
- precision
- structure solution

48